crea spine A C I F C A G E



INNOVATIVE

INNOVATIVE OLUTIONS

ACIF CAGE



The anterior cage implant is produced from PEEK material, a thermoplastic solution compatible with the body and causes no artifact. Completely compatible with MRI for in-vivo applications. The cage system is also available in forms of; rigid, bladed, expandable and bladed expandable alternatives:

- . Anatomically designed to match lordosis angle.
- Has a wide measurement range for height starting from 4 mm to 8 mm with 1 mm increment. External dimension range variations of 12x14mm and 14x16mm.
- . Has special threaded structure for extra stabilization in the application area.
- . The cage has large graft application space to provide maximum fusion.
- . The functional and robust application sets are designed and manufactures to facilitate the surgical application as well as to provide comfort for the patient.





ACIF CAGE INSTRUMENT SET

DESCRIPTION	QUANTITY	REF.NO
AWL	1	CS20001
CASPAR PIN DRIVER	1	CS20002
CASPAR PIN 12 mm	2	CS20003
CASPAR PIN 14 mm	2	CS20004
CASPAR PIN 16 mm	2	CS20005
CASPAR DISTRACTOR	1	CS20006
RASP 12x14	1	CS20007
RASP 14x16	1	CS20008
TRIAL HOLDER	1	CS20009
TRIAL 12x14 H4	1	CS20010
TRIAL 12x14 H5	1	CS20011
TRIAL 12x14 H6	1	CS20012
TRIAL 12x14 H7	1	CS20013
TRIAL 12x14 H8	1	CS20014
TRIAL 14x16 H4	1	CS20015
TRIAL 14x16 H5	1	CS20016
TRIAL 14x16 H6	1	CS20017
TRIAL 14x16 H7	1	CS20018
TRIAL 14x16 H8	1	CS20019
CERVICAL CAGE INSERTER AND REMOVER	1	CS20020
GREFT BOX	1	CS20021
GREFT PUNCH	1	CS20022
MALLET	1	CS20023
TRIAL AND IMPLANT BOX	1	CS20024
CONTAINER	1	CS20025
TRAY	1	CS20026

REF CODE	DESCRIPTION	ANGLE (S/I)
C01244P~C01248P		0/0 Degree
C01464P~C01468P		0/0 Degree
C01244L~C01248L		5/0 Degree
C01464L~C01468L	CERVICAL PEEK CAGE L 14x16 H: 4 ~ 8 mm	5/0 Degree

- Pin for X-ray and maximum hold
- Strong grip and easy remove
- Toothed inferior and superior surfaces
- Full compliance to posterior region
- Large graft space









